

2022 CATALOG ELISA Kits

Research Topics: Diabetes, Obesity, Metabolism, Cardiovascular Disease, Renal Disease, Liver Disease, Bone Metabolism, Inflammation

About Us

ImmunoDiagnostics Limited (IMD) is a spin-off biotech company from The University of Hong Kong, and has R&D activities in both Toronto, Canada and Hong Kong. The founders of IMD are academic professoriates with strong research background in immunology, metabolism, cardiovascular medicine, and antibody and protein engineering. The company has over 13 years of experience in biomarker discovery and development of highly-specific immunoassays for both research & *in vitro* diagnostics (IVD) for major infectious diseases, cardiometabolic disorders and autoimmune diseases.

IMD has established state-of-the-art platforms for expression, purification and functional characterization of bioactive proteins from different sources (E. coli, yeast, insect cells and mammalian cells), generation and validation of both polyclonal and monoclonal antibodies, identification and cloning of genes encoding monoclonal antibodies specific to a target of interest using Next Generation Sequencing. Furthermore, our immunoassays have been validated in a large number of unique clinical biobanks in Asia, Europe and North America. Thus far, IMD has developed several hundreds of research products, including bioactive proteins, validated antibodies and immunoassays, which have been widely used in over 30 countries, such as the United States, Canada, United Kingdom, China, Australia and Singapore. IMD products have been cited by many publications in prestigious journals including Cell, Cell Metabolism, Nature Series, JCI, PNAS.






Table of Contents


Insulin Related ELISA	4
#31380 Human Insulin ELISA Kit	5
#31780 Human C-peptide ELISA Kit	5
#33100 Wide Range Rat Insulin ELISA Kit	6
#32100 Wide Range Mouse Insulin ELISA Kit	6
#32393 High-sensitive Plus Mouse Insulin ELISA Kit	6
#33270 High-sensitive Rat Insulin ELISA Kit	6
#32270 High-sensitive Mouse Insulin ELISA Kit	6
#33380 Ultra-sensitive Rat Insulin ELISA Kit	6
#32380 Ultra-sensitive Mouse Insulin ELISA Kit	6
#36780 Mouse/Rat C-peptide ELISA Kit	7
#31420 Human Proinsulin ELISA Kit	8
Cytokines, Inflammation and Autoimmune Biomarker	9
#31220 High-sensitive Human C-reactive Protein ELISA Kit	9
#32750 Mouse IL-33 ELISA Kit	10
#32150 Mouse Major Urinary Protein-1 ELISA Kit	11
#32300 Mouse Proteinase 3 ELISA Kit	12
Bone Biomarker	13
#31109 Human P1NP ELISA Kit	13
Other Biomarker	14
#31990 Human PLAC9 ELISA Kit	14

Adipokines, Hepatokines & Other Hormones	15
#31010 Human Adiponectin ELISA Kit	15
#31012 Rapid Human Adiponectin ELISA Kit	16
#32010 Mouse Adiponectin ELISA Kit	17
#32180 Mouse FGF-21 ELISA Kit	17
#31180 Human FGF-21 ELISA Kit	18
#31200 Human FGF-19 ELISA Kit	19
#31770 Human Autotaxin ELISA Kit	20
#32770 Mouse Autotaxin ELISA Kit	21
#32700 Mouse PM20D1 ELISA Kit	21
#31700 Human PM20D1 ELISA Kit	22
#31390 Human Alpha-1 Antitrypsin ELISA Kit	23
#32980 Mouse GDF-15 ELISA Kit	24
#32280 Mouse CYSTM1 ELISA Kit	24
#31980 Human GDF-15 ELISA Kit	25
#31060 Human RBP4 ELISA Kit	26
#31030 Human FABP4 ELISA Kit	27
#31690 Human Galectin-3 ELISA Kit	28
#31400 Human Haptoglobin ELISA Kit	29
#31070 Human PAI-1 ELISA Kit	30
Liver Disease Biomarker	31
#31101 Human Thrombospondin-2 ELISA Kit	31
#31103 Human TIMP-1 ELISA Kit	32
Kidney Disease Biomarker	33
#31241 Rapid Human Cystatin C ELISA Kit	33
#32050 Mouse Lipocalin-2 ELISA Kit	34
#31050 Human Lipocalin-2 ELISA Kit	35



General Introduction to Insulin

Insulin is a polypeptide hormone predominantly produced by the pancreatic b-cells. It consists of A chain and B chain linked with one intra and two inter-disulfide bonds. Insulin production is initiated by the synthesis of a single-chain precursor called preproinsulin (*see diagram below*) in the b-cells of pancreatic islets. This precursor is translocated into the rough endoplasmic reticulum and the signal peptide of preproinsulin is cleaved to form proinsulin. In the final step, proinsulin is transported to golgi apparatus and C-peptide is proteolytically cleaved off to form insulin. The insulin formed has two polypeptide chains (A and B chain) with a molecular weight of approximately 5.8 kDa. It is responsible for maintaining the blood glucose levels. Measuring the concentration of insulin, proinsulin and C-peptide can help to assess bcell function as well as precision diagnosis of diabetes.



IMD has developed a comprehensive range of assays for quantitative measurement of insulin, C-peptide and proinsulin in different species with high sensitivity, high specificity and superior reproducibility. All these assays have been fully validated in preclinical and clinical samples. IMD insulin ELISA kits have a high customer demand, has been used by many academics and R & D researchers globally in the past decade.

IMD Insulin & C-peptide ELISA kits									
Cat. No. 31380	Human insulin ELISA	Cat. No. 31780	Human C-peptide ELISA						
Cat. No. 32270	High-sensitive mouse insulin ELISA	Cat.No. 33270	High-sensitive rat insulin ELISA						
Cat. No. 32393	High-sensitive plus mouse ELISA	Cat. No. 33100	Wide range rat insulin ELISA						
Cat. No. 32100	Wide range mouse insulin ELISA	Cat. No. 33380	Ultra-sensitive rat ELISA						
Cat. No. 32380	Ultra-sensitive mouse ELISA	Cat.No. 36780	Mouse/rat C-peptide ELISA						



Human Insulin & Human C-peptide ELISA Kits (Cat. No. 31380; 31780)

INTRODUCTION: Human insulin consists of two peptide chains, A chain (21 amino acids) and B-chain (30 amino acids) linked by disulfide bonds. The sequence of human insulin is as shown below:



C-peptide is a part of proinsulin which is proteolytically cleaved out during insulin formation. It consists of 31 amino acid residues (*see the sequence below*). The half -life of C-peptide (20-30 min) is longer than that of insulin (3-5 min). Therefore, C-peptide is a more stable parameter for assessing b -cell function.


Glu Ala Glu Asp Leu Gin Val Gly Gin Val Glu Leu Gly Gly Gly Pro Gly Ala Gly Ser Leu Gin Pro Leu Ala LeuGlu Gly Ser Leu Gin

IMD has developed monoclonal antibody-based ELISA for specific measurements of insulin and Cpeptide in human plasma, serum or other biological samples.

Assay performance & characteristics A. Typical representation of standard curve

Insulin standard curve (log-log)

Human C-peptide standard (4-parameter)



- Cheng, Sam Tsz Wai, et al. <u>PLoS One.</u>2016; 11(1):e0147391.
- Lin Z, Pan X, Wu F, et al. *Circulation*.2015; 131(21):1861-71.
- ✤ Wang, Lin, et al. <u>Pancreas.</u>2017; 46(3):395-404. ◆
- Li X, Cheng K. K, Liu Z, et al. <u>Nat Commun.</u>2016; 7:11740.
- Shu L, Hoo R L C, Wu X, et al. *Nat Commun.*2017; 8:14147.
- ✤ Wang, Lin, et al. <u>Pancreas.</u>2017; 46(3):395-404.



Mouse/Rat Insulin ELISA Kits (Cat. No. 32270; 32100; 32380; 32393; 33270; 33100; 33380)

INTRODUCTION: Mouse and rat produce two insulins, insulin I and insulin II due to the two-gene insulin system present in these rodents. Insulin consists of two peptide chains, A chain (21 amino acids) and Bchain (30 amino acids) linked by disulfide bonds. The amino acid sequence of insulins from both rodents are same however, two amino different in B chain of insulin I and II. The figure below represents mouse insulin 2B (Note: amino acids highlighted in red differ from human insulin).

	A-chai	n			5			s												
NH2-	Gly I	le Va	al Asp	Gin C	rs Cys	Thr	Ser	le Cys	Ser l	eu T	yr Gln	Leu G	lu Asn	Tyr Cys	Asn	соон				
	B-chai	n			s	P	ro Mo	use insulii	1B				s	S					Lys Mouse	insulin 1B
NH ₂	Phe V	al Ly	Gln	His Le	J Cys	Gly s	er His	Leu V	al Glu	Ala	Leu Ty	r Leu \	al Cys	Gly Glu	Arg Gly	Phe Phe	Tyr Ti	hr Pro	MetSer	соон

IMD has developed monoclonal antibody-based ELISA for specific measurements of insulin in mouse/rat plasma, serum or other biological samples.

Assay performance & characteristics (*high-sensitive mouse insulin*as an example***)**

Insulin (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.052	0
0.2	0.063	0.011
0.5	0.098	0.046
1.0	0.197	0.145
2.0	0.558	0.506
3.5	1.23	1.178
7	2.731	2.679

A. Typical representation of standard curve

B. Sensitivity: 0.2 ng/mL.


C. Precision: Intra-assay C.V < 6.1%

Inter-assay C.V. <2.9%

D. Recovery: >91%.

Cat. No.	Name of immunoassay kit	Assay range	Application
32393	High-sensitive plus mouse insulin	0.2-7 ng/mL	suitable for obese mice (sample volume: 5µL)
32100 33100	Wide range mouse insulin Wide range rat insulin	1 -50 ng/mL	suitable for cell culture studies
32380 33380	Ultra-sensitive mouse insulin Ultra-sensitive rat insulin	0.025-1.5 ng/mL	suitable for lean and healthy mice suitable for lean and healthy rat
33270	High-sensitive rat insulin	0.2-7 ng/mL	suitable for obese rat (sample volume :10µL)

*suitable for obese mice (sample volume :10µL)



Concentration (ng/mL)

6



Mouse/Rat C-peptide ELISA Kit (Cat. No. 36780)

INTRODUCTION: Mouse and rat produce two types of C-peptides, C-peptide I (C-I) and C-peptide II (C-II). Mouse C-I (29 amino acids) has three amino acids different from C-II (31 amino acids). However, both Cpeptides from rat consist of 31 amino acid residues with a difference in two amino acids (see the sequence below):

> Mouse C1: EVEDPQVEQLELGGSPG--DLQTLALEVARQ Mouse CII: EVEDPQVAQLELGGGPGAGDLQTLALEVAQQ


Rat C1: EVEDPQVPQLELGGGPEAGDLQTLALEVARQ Rat CII: EVEDPQVAQLELGGGPGAGDLQTLALEVARQ

IMD has developed monoclonal antibody-based ELISA for specific measurements of C-peptide I and II in mouse/rat plasma, serum or other biological samples.

C_peptide (ng/mL)	Absorbance (450 nm)	Blanked Absorbance	2.5 -	Mouse C-peptide standard curve
0	0.133	0	Ê 2.0-	/
0.3	0.151	0.018		/
0.6	0.203	0.07	(420 um) (420 um)	/
1.2	0.338	0.205		/
2.4	0.783	0.65	UB 1.0 -	/
6	2.479	2.346		
ensitivity : ().3 ng/ mļ		Absorbance	
	tra_assay C_V <4		01	1 10

Assay performance & characteristics A. Typical representation of standard curve (4-parameter)

С Inter assay C.V <7.5%.



References:

- ♦ So, W. Y., et al. *Cell Death Dis.*2015; 6:e1707.
- ◆ Liang Q, Zhong L, Zhang J, et al. *Diabetes*.2014; 63(12): 4064-4075.
- ✤ Wang, Lin, et al. <u>PLoS One.</u> 2015; 10(6):e0128216.
- Cheng K.K, Lam KS, et al. <u>Proc Natl Acad Sci U S A.</u> 2012; 109(23):8919-24.
- Cheng K.K, et al. <u>Diabetes.</u> 2014; 63(11): 3748-58.
- Lin Z, Tian H, et al. <u>Cell Metab.</u>2013; 17(5):779-89.
- Liang Q, Zhong L, Zhang J, et al. <u>Diabetes.</u>2014; 63(12): 4064-4075.
- Cheng, Sam Tsz Wai, et al. <u>PLoS One.</u>2016; 11(1):e0147391.

7



Human Proinsulin ELISA Kit

(Cat. No. 31420)

INTRODUCTION: Human proinsulin is a precursor of insulin which consists of A chain and B chain connected by C- peptide (*see diagram below*). It is a polypeptide composing of 86 amino acids, produced after cleavage of signal peptide from pre-proinsulin in pancreatic b-cells. Elevated proinsulin level in circulation is observed during aging process or in subjects with impaired insulin maturation. Additionally, patients with type 2 diabetes and metabolic syndrome also have higher circulating proinsulin levels than healthy individuals.

		A-chai	n																		
Gly	lle	Val	Glu	Gin Cy	Cys	Thr	Ser I	le C	ys Ser	Leu	Tyr G	In Leu	Glu	Asn Ty	r Cys	Asn			chain ceptide		
Arg																			chain		
Lys									C-pept	tide											
Gin	Leu S	ier Gh	y Glu	Leu A	la Leu	I Pro	Gin Le	u Ser	Gly A	la Gl	y Pro C	ily Gly	Gly	Leu Glu	Val	Gin (ily Val	Gin Leu	Asp G	lu Ala	Glu
																					Arg
		-					-									_	B-chai	n			Arg
Phe	Val	Asn G	iln H	is Leu	Cys G	ily Se	r His	Leu	Val Gl	u Ala	Leu T	yr Leu	Val 0	ys Gly	Glu	Arg G	ly Phe	Phe Tyr	Thr Pr	o Lys	Thr

IMD has developed monoclonal antibody-based ELISA for specific measurements of human proinsulin in plasma, serum or other biological samples.

Assay performance & characteristics A. Typical representation of standard curve

Human proinsulin (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.096	0
31.2	0.122	0.026
62.5	0.154	0.058
125	0.225	0.129
250	0.354	0.258
500	0.616	0.520
1000	1.214	1.118
2000	2.168	2.072

Diabetes serum sample	Concentration(pg/mL)
Sample 1	226
Sample 2	40

B. Sensitivity: 31.2 pg/ml

C. Precision:

Intra-assay Precision C.V.< 6.5%. Inter-assay Precision C.V.< 7.3%.

D. Recovery: 89.9%

E. Linearity: 98-113%



Human proinsulin standard (4-parameter)



High-sensitive Human CRP ELISA Kit (Cat. No. 31220)


INTRODUCTION: C-reactive protein (CRP) or Pentraxin 1, is a circulating protein mainly secreted by liver. It consists of five identical non-glycosylated subunits of 23 kDa, that give rise to a symmetrically arranged globular protein with molecular weight of approximately 120 kDa. In 2003, the Centre for Disease Control and Prevention (CDC) and the American Heart Association (AHA) issued a statement that identified CRP as the inflammatory marker best suited for use in current clinical practice to assess cardiovascular risk. Many epidemiologic studies have demonstrated that CRP is a strong independent predictor of future cardiovascular events, such as myocardial infarction, ischemic stroke, peripheral vascular disease, and sudden cardiac death. The CDC/ AHA guidelines support the use of CRP in primary prevention and set cut-off points according to relative risk categories: low risk (<1.0 mg/L), average risk (1.0-3.0 mg/L), and high risk (>3.0 mg/L).


IMD has developed a monoclonal antibody pair-based ELISA for quantitative determination of CRP in human plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve

CRP (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.066	0
0.312	0.125	0.059
0.625	0.183	0.117
1.25	0.317	0.251
2.5	0.545	0.479
5	0.955	0.889
10	1.665	1.599
20	2.498	2.432





B. Sensitivity: 0.312 ng/mL.

- C. Precision: Intra-assay C.V. 4.3%; Inter-assay C.V. 5.9%.
- D. Recovery: 98%.

E. Specificity: The antibody pair used in this assay is specific to human CRP and does not cross-react with mouse and rat CRP, and other cytokines or hormone molecules.

- 1. Chow, W.S. et al. (2013) J. Am. Heart Assoc. 2(1): e004176
- 2. Festa, A. *et al.* (2002) Diabetes **51**:1131-1137.
- 3. Verma, S. & Yeh, E.T. (2003) Am. J. Physiol. 285:R1253-R1258.
- 4. Jialal, I. et al. (2004) Hypertension 44: 6-11.
- 5. Liu, M.M. et al. (2009) Chest. 135(4):950-956




Mouse IL-33 ELISA Kit (Cat. No. 32750)

INTRODUCTION: Interleukin 33 (IL-33) is a nuclear cytokine composed of a precursor protein (166 amino acid) and hydrophobic signal peptide (26 amino acid). It is released from tissues/cells during damage or injury and binds to the receptors ST2 (also known as IL1RL1) and IL-1 Receptor Accessory Protein (IL1RAP), activating intracellular molecules in the NF-κB and MAP kinase signaling pathways that drive production of type 2 cytokines (e.g. IL-5 and IL-13) from polarized Th₂ cells. Altered IL-33 levels have been associated with several diseases such as asthma and allergy, and chronic inflammatory conditions.

IMD has developed monoclonal antibody-based ELISA for specific measurements of IL-33 in mouse plasma, serum or other biological samples.

Assay performance and characteristics



			2.5 -	
Mouse IL-	Absorbance	Blanked		
33 (pg/mL)	(450 nm)	Absorbance	2.0 -	7
0	0.11	0	Ē.	/
31.25	0.145	0.027	St 1.5 -	/
62.5	0.179	0.062	oe (/
125	0.236	0.118	- 0.1 Lpan	/
250	0.319	0.202	Absorbance (450 nm)	
500	0.65	0.532	₹ 0.5-	
1000	1.2	1.083	0.0	
2000	2.35	2.233		
•			10	100 1000
				Concentration (pg/mL)

A. Typical standard curve

B. Sensitivity: 6.5 pg/mL.


C. Precision: Intra-assay C.V. <6.5%.; Inter-assay C.V. <10%.


D. Spiking: The recovery rate of spiked samples ranges from 101.8-115.6%.

E. Linearity: The recovery ranges from 92.4-98.6%.

F. Validation:

Cell culture supernatants:





- 1. Schrader, J.W. (1986) Annu. Rev. Immunol. 4:205
- 2. Fung, M.C. et al. (1984) Nature 307:233
- 3. Zeyda, M. et al. (2013) Int. J. Obes. 37: 658-665
- 4. Miller, A.M. et al. (2010) Circ. Res. 107(5): 650-8
- 5. Giannoudaki, E. *et al.* (2019) Nat. Commun. **10**:4003
- 6. Duffen, J. et al. (2018) J. Immunol. 200(4):1347-1359




Mouse Major Urinary Protein-1 (MUP1) ELISA Kit (Cat. No. 32150)

INTRODUCTION: Major urinary protein 1 (Mup1), also known as Mup7, Up-1, Ltn-1, Mup-1, Mup-a, Mup10 and Lvtn-1, is a low molecular weight secretory protein produced predominantly from the liver. Structurally it belongs to the lipocalin family, which carries small hydrophobic ligands such as pheromones. It has been demonstrated that Mup1 is an important player in regulating energy expenditure and metabolism in mice. Circulating Mup1 is markedly decreased in db/db mice with obesity and diabetes. MUP1 treatment decreases glucose intolerance as well as insulin resistance in obese/diabetic mice.

30-

Assay performance & characteristics



	5.0 7			
	د	Blanked Absorbance	Absorbance (450 nm)	Mup1 (ng/mL)
	LL 20-	0	0.071	0
	(420 nm)	0.046	0.117	0.625
/	Absorbance	0.103	0.174	1.25
_	-0.1 LDa	0.208	0.279	2.5
	osqv 0.5	0.441	0.512	5
	▲]	0.852	0.927	10
	0.0 -	1.579	1.65	20
1 10	0.1	2.484	2.555	40
Concentration (ng/mL)	1	•	•	

Typical standard curve

Sensitivity: 0.2 ng/mL.; Intra-assay C.V. <5.3%.; Inter-assay C.V. <4.4%

References for MUP-1:

- 1. McIntosh, I. and Bishop, J.O. (1989) Mol. Cell. Biol. 9:2202-2207.
- 2. Hui, X. et al. (2009) J. Biol. Chem. 284(21):14050-7.
- 3. Zhou, Y. et al. (2009) J. Biol. Chem. 284(17):11152-9.
- 4. Chen, C.C. et al. (2015) Biochem. Biophys. Res. Commun. 460(4):1063-8.



Mouse Proteinase 3 (PR3) ELISA kits (Cat. No. 32300)

INTRODUCTION: Proteinase 3 (PR3), also known as myeloblastin, Wegener autoantigen, azurophil granule protein-7 or p29b, is one of the hematopoietic serine proteases localized in the primary granules of polymorphonuclear neutrophils (PMNs). It is a 29 kDa glycoprotein made of 222 amino acids and is released during neutrophilic inflammation. The primary function of PR3 is to participate in direct intracellular killing of phagocytosed pathogens in phagolysosomes and degradation of extracellular matrix components at inflammatory sites. PR3 has also been shown to process some pro-inflammatory cytokines, activate mitogen activated protein kinase (MAPK) signalling through proteinase activated receptor-1 (PAR1), and induce endothelial cell apoptosis through NF-κB signalling pathways. PR3 is identified as the target autoantigen of anti-neutrophil cytoplasmic autoantibodies (ANCA) in Wegener granulomatosis. Increased PR3 levels have been reported in patients with acute myocardial infarction, and in subjects with type 1 diabetes.

Mouse PR3 standard curve (4-parameter)

Assay performance & characteristics

	1		
Mouse PR3	Absorbance	Blanked	25-
(ng/mL)	(450 nm)	Absorbance	(u 20 - 0
0	0.08	0	(420
0.1	0.117	0.037	Approximation of the second se
0.2	0.154	0.074	20 1.0 -
0.4	0.243	0.163	• •
0.8	0.416	0.336	
1.6	0.766	0.686	0.0 -
3.2	1.483	1.403	0.1
6.4	2.604	2.524	Concentration (ng/mL)

Typical standard curve

Sensitivity: 0.1 ng/mL.; Intra-assay C.V. <10%.; Inter-assay C.V. <10%



Human P1NP ELISA Kit (Cat. No. 31109)

INTRODUCTION: Procollagen type 1 N-terminal propeptide (P1NP) is the N-terminal extension of procollagen type 1 found in bone. During collagen formation, this N-terminal is proteolytically removed along with the C-terminal extension. P1NP is usually secreted into blood stream as unstable trimeric structure which subsequently undergo thermal degradation to form a stable monomeric structure. The International Osteoporosis Foundation recommends P1NP as a bone formation biomarker. Various studies of hormone replacement therapy in postmenopausal women portrayed the significance of P1NP as biomarker for bone formation. Additionally, P1NP was observed as good indicator for fibrogenesis in rats and osteoblastic metastases in prostate carcinoma patients.


IMD has developed a monoclonal antibody pair-based ELISA for specific measurements of P1NP in human plasma, serum or other biological samples.

Assay performance & characteristics

Human P1NP (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.094	0
62.5	0.133	0.039
125	0.187	0.093
250	0.293	0.199
500	0.474	0.38
1000	0.812	0.718
2000	1.317	1.223
4000	2.096	2.002

A. Typical standard curve

Human P1NP standard curve (4-parameter)



B. Sensitivity: 62.5 pg/mL.

C. Precision: Intra-assay C.V.< 3.3%; Inter-assay C.V. < 5.1%

- 1. Garnero, P. et al. (2008) Clin. Chem. 54(1):188-196
- 2. Kuo, T.R. and Chen, C.H. (2017) Biomark. Res. 5:18
- 3. Hannon, R. et al. (2009) J. Bone Miner. Res. 13:7
- 4. Eastell, R. et al. (2005) Curr. Med. Res. Opin. 22(1):61-6
- 5. Vasikaran, S. et al. (2010) Osteoporos Int. 22(2):391-420
- 6. Funck-Brentano, T. et al. (2011) Semin. Arthritis Rheum. 41(2):157-69




Human PLAC9 ELISA Kit (Cat. No. 31990)

INTRODUCTION: Placenta-specific protein 9 (PLAC9) is a putative secreted protein highly enriched in placenta. It plays a major role in embryonic development. During human embryogenesis, PLAC9 levels were upregulated at 8-9 weeks. PLAC9 was found to be involved in protein interactions in human liver. Overexpression of PLAC9 inhibited cell proliferation capacity in hepatic cells by modifying cell cycle related proteins.

Assay performance & characteristics

A. Typical standard curve

Human PLAC9 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.108	0
0.156	0.126	0.018
0.312	0.157	0.049
0.625	0.306	0.198
1.25	0.648	0.540
2.5	1.107	0.999
5	1.584	1.476
10	2.192	2.084



Human PLAC9 standard curve

B. Sensitivity: 0.156 ng/mL.

- 1. Galaviz-Hernandez C., et al., (2003) Gene, 309(2):81-9
- 2. Yi H., et al., (2010) FASEB J., 24(9): 3341–3350
- 3. Xue L., et al., (2011) Int J Biol Sci, 7(7):1068-1076
- 4. Wang J., et al., (2011) Mol Syst Biol., 7: 536



Human Adiponectin ELISA Kit (Cat. No. 31010)


INTRODUCTION: Adiponectin, also known as apM1, Acrp30, adipoQ and GBP28, is a circulating hormone predominantly produced from adipose tissue. It is a glycoprotein with molecular weight of approximately 30 kDa and have circulating levels of $2-30 \,\mu$ g/mL in humans. Structurally it resembles complementary factor C1q. During assembly, adiponectin form complexes of different molecular weight and function. Adiponectin has anti-diabetic, insulin-sensitizing and anti-inflammatory properties. Decreased circulating levels of plasma adiponectin (hypoadiponectinemia) are associated with increased body mass index (BMI), and decreased insulin sensitivity. A large number of longitudinal studies in different ethnic groups uniformly demonstrated that circulating adiponectin levels are decreased significantly in type 2 diabetes (T2D) and related complications, and low circulating adiponectin levels indicate the increased risk for the development of T2D.

IMD has developed monoclonal antibody-based ultrasensitive ELISA for specific measurements of adiponectin in human plasma, serum or other biological samples.

Assay performance & characteristics

Adiponectin	Absorbance	Blanked
(ng/mL)	(450 nm)	Absorbance
0	0.068	0
1.56	0.112	0.044
3.12	0.156	0.088
6.25	0.248	0.18
12.5	0.418	0.35
25	0.724	0.656
50	1.384	1.316
100	2.469	2.401

A. Typical standard curve



B. Sensitivity: 1.56 ng/mL.

C. Precision: Intra-assay C.V.< 4.18%.; Inter-assay C.V.< 5.01%.

D. Recovery: 97-105%.

E. Specificity:

The antibody pair used in this assay is specific to human adiponectin and does not cross-react with mouse and rat adiponectin, and other cytokines or hormone molecules including human resistin, $TNF\alpha$, ANGPTL4, insulin, leptin and IL6.

- 1. Trujillo, M.E. and Scherer, P.E. (2005) J. Intern. Med. 257:167-175
- 2. Robinson, K. *et al.* (2011) Crit. Care **15(2**):221
- 3. Berg, A.H. et al. (2001) Nat. Med. 7:947-953.
- 4. Yamauchi, T. et al. (2001) Nat. Med. 7:941-946.
- 5. Xu, A. *et al.* (2003) J. Clin. Invest. **112**:91-100.
- 6. Lin, Z. et al. (2013) Cell. Metab. 17(5):779-89


Rapid Human Adiponectin ELISA Kit (Cat. No. 31012)

INTRODUCTION: Adiponectin, also known as apM1, Acrp30, adipoQ and GBP28, is a circulating hormone predominantly produced from adipose tissue. It is a glycoprotein with molecular weight of approximately 30 kDa and have circulating levels of 2-30 µg/mL in humans. Structurally it resembles complementary factor C1q. During assembly, adiponectin form complexes of different molecular weight and function. Adiponectin has anti-diabetic, insulinsensitizing and anti-inflammatory properties. Decreased circulating levels of plasma adiponectin (hypoadiponectinemia) are associated with increased body mass index (BMI), and decreased insulin sensitivity. A large number of longitudinal studies in different ethnic groups uniformly demonstrated that circulating adiponectin levels are decreased significantly in type 2 diabetes (T2D) and related complications, and low circulating adiponectin levels indicate the increased risk for the development of T2D.

IMD has developed monoclonal antibody-based ELISA for rapid, specific measurements of adiponectin in human plasma, serum or other biological samples within 90 minutes.

Assay performance & characteristics A. Typical standard curve

Adiponectin (ng/mL)	Absorbance (450 nm)	Blanked absorbance
0	0.056	0
7.8	0.105	0.049
15.6	0.158	0.102
31.2	0.255	0.199
62.5	0.463	0.407
125	0.845	0.789
250	1.608	1.552
500	2.599	2.543



B. Sensitivity: 7.8 ng/mL.

C. Precision: Intra-assay C.V. < 4.7%; Inter-assay C.V. < 6.9%

D. Recovery: 91%.

E. Specificity:


The antibody pair used in this assay is specific to human adiponectin and does not cross-react with mouse and rat adiponectin, and other cytokine or hormone molecules including human resistin, $TNF\alpha$, ANGPTL4, insulin, leptin and IL6.

- 1. Trujillo, M.E. and Scherer, P.E. (2005) J. Intern. Med. 257:167-175
- 2. Robinson, K. *et al.* (2011) Crit. Care **15(2**):221
- 3. Berg, A.H. et al. (2001) Nat. Med. 7:947-953.
- 4. Yamauchi, T. et al. (2001) Nat. Med. 7:941-946.
- 5. Xu, A. et al. (2003) J. Clin. Invest. 112:91-100.
- 6. Lin, Z. et al. (2013) Cell. Metab. 17(5):779-89

ImmunoDiagnostics

Mouse Adiponectin ELISA Kit (Cat. No. 32010)

INTRODUCTION: Refer Cat. No. 31010 for more information. **Assay performance & characteristics Typical standard curve**



Sensitivity: 2 ng/mL; Intra-assay C.V. <7.6%., Inter-assay C.V. <6.4%.; Recovery: 95%


References for Adiponection Assay:

- 1. Xu, A. et al. (2004) Endocrinology 145(2):487-94.
- 2. Xu, A. et al. (2005) J. Biol. Chem. 280(18): 18073-80.
- 3. Palanivel, R. et al. (2007) Cardiovasc. Res. 75(1): 148-57.
- 4. Wong, J.T. et al. (2011) Cell Metab. 14(1):104-15.
- 5. Lin, Z. et al. (2014) J. Hepatol. 61(4):825-31.
- 6. Hui, X. et al. (2015) Cell Metab. 22(2):279-90.

Mouse/Rat FGF-21 ELISA Kit (Cat. No. 32180)

INTRODUCTION: Refer Cat. No. 31180 for more information. **Assay performance & characteristics**

Mouse FGF-21 standard curve (4-parameter)



Sensitivity: 31.2 pg/mL; Intra-assay C.V. <5.6%, Inter-assay C.V. <7.1%.; Recovery: 96%.

References for FGF-21 Assay:

- 1. Chen, W. et al. (2011) J. Biol. Chem. 286(40):34559-66.
- 2. Li, H. et al. (2012) Diabetes 61(4):797-806.
- 3. Lin, Z. et al. (2013) Cell Metab. 17(5):779-89.
- 4. Huang, Z. et al. (2017) Cell Metab. 26(3):493-508.
- 5. Pan, X. et al. (2018) Cell Metab. 27(6):1323-1337

Website: www.immunodiagnostics.com.hk E-mail: sales@immunodiagnostics.com.hk / info@immunodiagnostics.ca Tel: +852 3502 2780/+1-437-886-5136



Human FGF-21 ELISA Kit (Cat. No. 31180)

INTRODUCTION: Fibroblast growth factor 21 (FGF-21) is an endocrine member of the FGF superfamily produced predominantly from the liver. Additionally, it is also expressed in several other metabolically active organs such as adipose tissues, pancreatic islets and skeletal muscle. FGF-21 is a potent metabolic regulator with therapeutic benefits on obesity and its related metabolic complications, partly by induction of adiponectin. Paradoxically, circulating levels of FGF-21 are elevated in obesity, type-2 diabetes and related complications, possibly due to FGF-21 resistance or compensatory responses. Prospective studies in different ethnic groups have identified high level of serum FGF-21 as an independent predictor for type-2 diabetes and atherosclerosis.


IMD has developed a highly sensitive ELISA for specific measurements of FGF-21 in human plasma, serum or other biological samples, which have been used widely for clinical research.

Assay performance & characteristics

A. Typical standard curve

Human FGF-21 (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.094	0
30	0.127	0.033
60	0.158	0.064
120	0.224	0.13
240	0.362	0.268
480	0.626	0.532
960	1.141	1.047
1920	2.14	2.046

Human FGF-21 standard curve (4-parameter)



B. Sensitivity: 30 pg/mL.

D. Specificity:

The recombinant protein does not cross-react with mouse FGF-21, human FABP4, human LCN2, human adiponectin and human ANGPL4

E. Linearity:

	Measured	Expected	Recovery
Dilution	(pg/mL)	(pg/mL)	(%)
None	448	448	100
1/2	245	224	109
1/4	123	112	109
1/8	63.6	56	113

References:

1. Zhang, X. (2008) Diabetes 57(5):1246-53.

- 2. Yu, H. (2011) Clin. Chem. 57(5):691-700
- 3. Li, H. J. (2013) Hepatol. 58(3):557-63
- 4. Ong, K.L. (2015) Diabetologia. 58(9):2035-44.
- 5. Lee, C.H. (2017) J. Am. Heart Assoc. 6(6): e005344

C. Precision: Intra-assay C.V. < 5.1%; Inter-assay C.V. < 4.6%



Human FGF-19 ELISA Kit (Cat. No. 31200)

INTRODUCTION: FGF-19 and its mouse ortholog FGF-15 is a hormone-like molecule predominantly expressed in the distal small intestine. It is an important regulator of bile acid homeostasis, through its actions in the liver, where it binds to its receptor complex comprising of FGFR4 and the single transmembrane co-receptor β -klotho. FGF-19 has also been shown to possess potent beneficial effects on glucose metabolism. Transgenic mice expressing FGF-19 had increased energy expenditure and improved glucose tolerance, and administration of recombinant FGF-19 protein prevented the development of glucose intolerance in both high-fat-fed mice and leptin-deficient mice. Decreased circulating level of FGF-19 is observed in patients with metabolic syndrome, impaired glucose tolerance and type-2 diabetes.


IMD has developed a highly specific antibody-based ELISA for accurate measurements of FGF-19 in human plasma, serum or other biological samples.

Assay performance & characteristics

Human FGF-19 (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.088	0
31.2	0.121	0.033
62.5	0.155	0.067
125	0.225	0.137
250	0.374	0.286
500	0.66	0.572
1000	1.201	1.113
2000	2.136	2.048

A. Typical standard curve

Human FGF-19 standard (4-parameter)



B. Sensitivity : 31.2 pg/mL.

C. Precision: Intra-assay C.V.< 4.5%; Inter-assay C.V.< 5.6%.

D. Specificity:

The antibodies used in this assay are specific to human FGF-19 and do not cross-react with human adiponectin, FGF-21, FABP4, LCN2, RBP4 and PAI-1.

- 1. Dolegowska, K. et al. (2019) J. Physiol. Biochem. 7(5):229-240
- 2. Nishimura, T. et al. (1999) Biochem. Biophys. Acta. 1444:148
- 3. Degirolamo, C. et al. (2015) Nat. Rev. Drug Discov. 15:51-69
- 4. Beenken, A. et al. (2009) Nat. Rev. Drug Discov. 8: 235-253.
- 5. Zhang, J. et al. (2017) Sci. Rep. 7(1):6042.
- 6. Fang, Q. et al. (2013) Diabetes Care 36(9): 2810-2814.




Human Autotaxin ELISA Kit (Cat. No. 31770)

INTRODUCTION: Autotaxin, also known as ENPP2 (ectonucleotide pyrophosphatase/phosphodiesterase 2), is a secreted glycoprotein which belongs to the ectonucleotide pyrophosphatase/phosphodiesterase (NPP) family. Autotaxin exhibits unique lysophospholipase D activity which allows formation of lysophosphatidic acids (LPA) and choline from lysophosphatidylcholine. Autotaxin is highly expressed in adipose tissues and its expression is markedly elevated during adipogenesis. Increased circulating level of autotaxin is associated with obesity-related metabolic complications. Additionally, elevated autotaxin levels are observed in patients with breast and lung carcinomas, tumor metastasis and liver fibrosis.

IMD has developed a highly specific-ELISA for accurate measurements of autotaxin in human plasma, serum or other biological samples.

Assay performance & characteristics

Human autotaxin (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.061	0
1.56	0.141	0.08
3.12	0.222	0.161
6.25	0.368	0.307
12.5	0.656	0.595
25	1.191	1.13
50	2.355	2.294



A. Typical standard curve

B. Sensitivity: 0.78 ng/mL.

C. Precision: Intra-assay C.V. <7%.; Inter-assay C.V. <7%.

D. Spiking: The recovery of human autotaxin spiked to 0, 5, 10 ng/mL falls in the range 90-110%.

E. Linearity:

Dilution	Measured ng/mL	Expected ng/mL	Recovery %
	1284	1284	100
1:40	856.8	856.8	100
	1151	1151	100
	1337.1	1284	104.1
1:80	881	856.8	102.8
	1216	1151	105.6

- 1. Cimpean, A. et al. (2017) Biochem. J. 474(24): 4273.
- 2. Okudaira, S. (2010) Biochimie 92(6): 698-706.
- 3. Umezu-Goto, M. et al. (2002): J. cell biol. 158(2): 227-233.
- 4. D'Souza, K. et al. (2018) J. Lipid Res. 59(10):1805-1817.
- 5. Lee, D. et al. (2018) Cancer Metastasis Rev. 31(2-3):509-518.
- 6. Brandon, J.A. et al. (2019) PLoS One 7:14(2).




Mouse Autotaxin ELISA Kit (Cat. No. 32770)

INTRODUCTION: Refer Cat. No. 31770 for more information. & 31700

Assay performance & characteristics

Typical standard curve



Sensitivity:0.312 ng/mL; Intra-assay C.V. <10%, Inter-assay C.V. <10%; Recovery: 80-120%.

Mouse PM20D1 ELISA Kit (Cat. No. 32700)

INTRODUCTION: Refer Cat. No. 31770 for more information.


Assay performance & characteristics

Typical standard curve

Mouse PM20D1 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.075	0
0.312	0.107	0.032
0.625	0.141	0.066
1.25	0.205	0.13
2.5	0.35	0.275
5	0.596	0.521
10	1.189	1.114
20	2.054	1.979

Sensitivity: 0.156 ng/mL; Intra-assay C.V. <3.8% Inter-assay C.V. <4.9%; Recovery: 93%-115%

Mouse PM20D1 standard curve (4-parameters)





Human PM20D1 ELISA Kit (Cat. No. 31700)

INTRODUCTION: Peptidase M20 domain containing 1 (PM20D1), also known as bidirectional N-fatty-acyl amino acid synthase/hydrolase, regulates the production of N-fatty-acyl amino acids. In brown and beige adipocytes, these metabolites may act as endogenous chemical uncouplers of mitochondrial respiration in an uncoupling protein 1-independent manner. Mice with increased circulating PM20D1 have augmented respiration and increased N-acyl amino acids in blood. Lastly, administration of N-acyl amino acids to mice improves glucose homeostasis and increases energy expenditure. Natural human genetic variation determines basal and inducible expression of PM20D1, an obesity-associated gene. PM20D1 is a quantitative trait locus associated with Alzheimer's disease.

IMD has developed a highly specific antibody-based ELISA for quantitative measurements of PM20D1 in human plasma, serum or other biological samples.

Assay performance & characteristics


Human PM20D1 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.105	0
0.39	0.131	0.03
0.78	0.139	0.043
1.56	0.162	0.079
3.125	0.25	0.202
6.25	0.505	0.404
12.5	0.846	0.809
25	1.74	1.615

A. Typical standard curve

B. Sensitivity: 0.209 ng/mL. C. Precision: Intra-assay C.V. < 3.7%; Inter-assay C.V. < 3.8% D. Spike:

Spike level	Expected (ng/mL)	Observed (ng/mL)	Recovery (%)
Low spike (1.25 ng/mL)	1.25	1.22	98
Medium spike (5 ng/mL)	5	5.32	106.5
High spike (20 ng/mL)	20	18.65	93.2

- 1. Long, J.Z. *et al.* (2016) Cell **166(2)**:424-435
- 2. Sanchez-Mut, J.V. (2018) Nat. Med. **24(5):**598-603
- 3. Long, J.Z. et al. (2018) PNAS **115(29**): E6937-E6945
- 4. Benson, K.K. *et al.* (2019) PNAS **116(46**):23232 23242



E. Specificity : No cross reactivity with recombinant
mouse PM20D1 protein.
F. Linearity:

Dilution	Measured (ng/ml)	Expected (ng/ml)	Recovery (%)
1/2	10.11	9.51	106.3
1/4	4.91	4.9	100.2
1/8	2.56	2.69	95.2



Human Alpha-1 Antitrypsin (A1AT) ELISA Kit (Cat. No. 31390)

INTRODUCTION: Alpha-1 antitrypsin (A1AT), also known as alpha-1 protease inhibitor or SERPINA1, is a 52 kDa serine protease inhibitor, which is mainly synthesized in the liver and is also produced in monocytes, macrophages, dendritic cells, etc. The primary function of A1AT is to inhibit the actions of proteolytic enzymes, such as neutrophil elastase (NE), proteinase 3 (PR3) and cathepsin G (CG), and to provide essential protection to host tissues from non-specific injury during periods of inflammation. A1AT has also been reported to regulate neutrophil chemotaxis, enhance insulin secretion and protect β -cells against cytokine-induced apoptosis. It possesses anti-inflammatory as well as immunomodulatory properties. Systemic deficiency of A1AT due to genetic mutations can result in a number of diseases, such as chronic obstructive pulmonary disease (COPD), systemic sclerosis, liver injury, cirrhosis and hepatocellular carcinoma. Furthermore, A1AT disorder has also been known to be involved in the development of diabetes mellitus.


IMD has developed a highly specific antibody-based ELISA for specific measurements of A1AT in human plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve

Human A1AT (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.08	0
0.312	0.129	0.049
0.625	0.166	0.086
1.25	0.267	0.187
2.5	0.445	0.365
5	0.816	0.736
10	1.47	1.39
20	2.55	2.47

Human A1AT standard curve (4-parameter)



B. Sensitivity: 0.312 ng/mL

C. Precision: Intra-assay C.V. 4.9%, inter-assay C.V. 5.3%

- 1. Janciauskiene, S.M. et al. (2011) Respir. Med. 105:1129-1139.
- 2. Miravitlles, M. (2012) Curr. Opin. Pharmacol. 12:309-314.
- 3. Wang, Y. et al. (2014) Diabetes 63(12):4239-48.
- 4. Sandstrom, C.S. et al. (2008) Diabet. Med. 25(11):1370-3
- 5. Weir, G.C. et al. (2018) Pediatr. Diabetes 19(5):945-954
- 6. Mitchell, E.L. & Khan, Z. (2018) Curr. Pathobiol. Rep. 6(1):97
- 7. Takei, N. et al. (2019) Int. J. Chron. Obstruct. Pulmon. Dis. 14:2885-2893




Mouse GDF-15 ELISA Kit (Cat. No. 32980)

INTRODUCTION: Refer Cat. No. 31980 for more information about GDF-15.

Assay performance & characteristics Typical standard curve

Mouse GDF-15 (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.091	0
3.9	0.126	0.035
7.8	0.164	0.073
15.6	0.246	0.155
31.2	0.404	0.313
62.5	0.692	0.601
125	1.259	1.168
250	2.216	2.125



Sensitivity: 3.9 pg/mL; Intra-assay C.V. < 2.9%, Inter-assay C.V. < 6.1%.; Recovery: 82-101%

Mouse CYSTM1 ELISA Kit (Cat. No. 32280)


INTRODUCTION: Cysteine-rich transmembrane module containing 1 (CYSTM1) also known as C5orf32, ORF1-FL49, belongs to CYSTM family which is a part of tail-anchored membrane proteins in eukaryotes. It consists of 104 and 97 amino acids in mouse and human, respectively. CYSTM1 is a transmembrane protein and based on Gene Ontology (GO) the function of CYSTM1 is classified as neutrophil degranulation. CYSTM1 was also reported as one of the candidate biomarkers for Huntington's disease.

Assay performance & characteristics

Typical standard curve

CYSTM1	Absorbance	Blanked
(ng/mL)	(450 nm)	Absorbance
0	0.137	0
0.156	0.162	0.025
0.312	0.194	0.057
0.625	0.273	0.136
1.25	0.429	0.292
2.5	0.71	0.573
5	1.316	1.179
10	2.458	2.321

Mouse CYSTM1 standard curve (4-parameter)



Sensitivity: 0.156 ng/mL; Intra-assay C.V. < 5%, Inter-assay C.V. <5.5%.; Recovery: 82-101%

References:

1. Mastrokolias, A, et al. (2015) Eur. J. of Hum Genet. 23(10): 1349-1356.

2. Xu, Y. et al. (2018) Plant Cell Physiol. 59(2): 423-438.

Mouse GDF-15 standard curve (4-parameter)



Human GDF-15 ELISA Kit (Cat. No. 31980)


INTRODUCTION: Growth differentiation factor 15 (GDF-15), also known as macrophage inhibiting cytokine 1 (MIC-1), placental transformation growth factor (PTGF- β), prostate derived factor (PDF), placental bone morphogenetic protein (PLAB), NSAID activated gene-1 (NAG-1) and PL74, belongs to the transforming growth factor β superfamily. GDF-15 is synthesized as a 62 kDa precursor protein, which, after cleavage by furin-like protease, is secreted as 25-kDa disulfide-linked dimer. GDF-15 is an important regulator of appetite and energy homeostasis. It exerts its effects via its receptor called glial-derived neurotrophic factor (GDNF) receptor alpha-like (GFRAL). GDF-15 expression and serum levels increase in response to many stimuli that initiate cell stress and as a part of a wide variety of disease processes, most prominently cancer, diabetes and cardiovascular disease. Circulating GDF-15 is a promising prognostic marker for heart failure and cancer.

IMD has developed highly sensitive monoclonal antibody-based ELISA for specific measurements of GDF-15 in human plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve

Human GDF-15 (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.074	0
7.8	0.13	0.056
15.6	0.197	0.123
31.2	0.338	0.264
62.5	0.577	0.503
125	1.017	0.943
250	1.719	1.645
500	2.626	2.552



B. Sensitivity: 7.8 pg/mL.

C. Precision: Intra-assay C.V. <10%.; Inter-assay C.V. <10%.

- 1. Adela, R. & Banerjee, S.K. (2015) J. Diabetes Res. 2015:1-14
- 2. Baek, S.J. & Eling, T. (2019) Pharmacol. Ther. 198:46-58
- 3. Bootcov, M.R. et al. (1997) Proc. Natl. Acad. Sci. U S A 94:11514-11519.
- 4. Welsh, J.B. *et al.* (2003) Proc. Natl. Acad. Sci. U S A **100**:3410-3415.
- 5. Chrysovergis, K. et al. (2014) Int. J. Obes. 38:1555-1564
- 6. Wang, T.J. *et al.* (2012) Circulation **126**:1596-1604.
- 7. Emmerson, P.J. et al. (2017) Nat. Med. 23:1215-1219
- 8. Nair, V. et al. (2017) J. Am. Soc. Nephol. 28:2223-2240.
- 9. Hsu, J.Y. et al. (2017) Nature 550: 255-259.




Human RBP4 ELISA Kit (Cat. No. 31060)

INTRODUCTION: Retinol binding protein 4 (RBP4) is a member of lipocalin family which transport small hydrophobic molecules. It is an adipocyte-secreted adipokine functionally involved in the transport of retinol and insulin resistance. RBP4 has a molecular mass of ~21kDa and composed of 201 amino acids. Half-life of RBP4 in serum is increased due to the interaction of RBP4-retinol complex with prealbumin and transthyretin. Elevated circulating level of RBP4 is observed in obese individuals and patients with insulin resistance, type 2 diabetes and atherosclerosis.

IMD has developed monoclonal antibody pair-based ELISA for specific measurements of RBP4 in human

plasma, serum or other biological samples.

Assay performance & characteristics



A. Typical standard curve

B. Sensitivity: 3.12 ng/mL.

C. Precision: Intra-assay C.V. < 3.63%; Inter-assay C.V. < 4.35%.

D. Recovery: 92-115%

E. Specificity:

The antibody pair used in this assay is specific to human RBP4 and does not cross-react with mouse and rat RBP4, and other cytokines or hormone molecules.

- 1. Ong, K.L. (2012) J. Clin. Endocrinol. Metab. 97(12):4701-8.
- 2. Zabetian-Targhi, F. et al. (2015) Adv. Nutr. 6(6):748-762
- 3. Xiao, Y. (2013) PLoS One 8(6): e66607.
- 4. Graham, T.E. et al. (2006) N. Engl. J. Med. 354: 2552–2563
- 5. Ingelsson, E. et al. (2009) Diabetes Care 32: 733-735
- 6. Wong, Y.K. *et al.* (2018) Arterioscler Thromb Vasc Biol. **38(10)**:2519-2527.
- 7. Jüllig, M. *et al.* (2014) PLoS One **9(5)**: e96489.
- 8. Alkharfy, K.M. et al. (2012) PLoS One 7(10): e48612



Human FABP4 ELISA Kit (Cat. No. 31030)


INTRODUCTION: Fatty-acid binding protein 4 (FABP4), also known as adipocyte fatty-acid binding protein (A-FABP) or adipocyte protein 2 (aP2), is a member of intracellular lipid chaperone. It is highly expressed in adipocytes and has a molecular weight of approximately 15 kDa. FABP4 is secreted from adipose tissue into bloodstream and act as a pro-inflammatory adipokine contributing to obesity-related cardiometabolic syndrome. Serum level of FABP4 increases during obesity, metabolic syndrome, non-alcoholic fatty liver disease and other cardiovascular disease. Increased expression of FABP4 is also observed during the progression of various cancers. Furthermore, high serum FABP4 is a predictor for type 2 diabetes and cardiovascular mortality.

IMD discovered the circulating form of FABP4 and has developed a monoclonal antibody-based ELISA for specific measurements of FABP4 in human plasma, serum or other biological samples.

Assay performance & characteristics

Absorbance	Blanked
(450 nm)	Absorbance
0.1	0
0.158	0.058
0.228	0.118
0.332	0.232
0.554	0.454
0.957	0.857
1.547	1.447
2.289	2.189
	(450 nm) 0.1 0.158 0.228 0.332 0.554 0.957 1.547

A. Typical standard curve



B. Sensitivity: 0.39 ng/mL.

C. Precision: Intra-assay C.V.< 4.1%; Inter-assay C.V.< 4.5%

D. Specificity:

The antibodies used in this assay are specific to human FABP4 and do not cross-react with mouse and rat FABP4, and other cytokines or hormone molecules.

- 1. Furuhashi, M. et al. (2014) Clin. Med. Insights Cardiol. 8:23-33
- 2. Lee, C.H. et al. (2018) Clin. Chem. 64(10):1496-1504.
- 3. Xu, A. et al. (2006) Clin. Chem. 52(3):405-13.
- 4. Xu, A. et al. (2007) Circulation 115:1537-1543.
- 5. Rhee, E.J. *et al.* (2009) Eur. J. Endocrinol. **160(2)**:165-72.
- 6. Tso, A.W. *et al.* (2007) Diabetes Care **30(10)**:2667-72
- 7. Hyun, K. J. et al. (2009) Diabetes Care **32(1)**: 147-152.



Human Galectin-3 ELISA Kit (Cat. No. 31690)


INTRODUCTION: Galectin-3 is a member of the β -galactoside-binding lectins family which has conserved carbohydrate recognition domain (CDR). It is a 35 kDa protein located in cytoplasm, nucleus and extracellular spaces. It is also also a member of the beta-galactoside-binding protein family that plays an important role in cell-cell adhesion, cell matrix interactions, macrophage activation, angiogenesis, apoptosis and insulin resistance. Elevated serum galectin-3 levels were observed in Behcet's disease, thyroid, Alzheimer's disease, cardiovascular disease such as left atrial appendage (LAA) stroke and in several types of cancers especially when it is metastatic. Moreover, circulating levels of galectin-3 were higher in obese patients and it is an indicator of insulin resistance.

IMD has developed a highly specific antibody-based ELISA for quantitative measurements of galectin-3 in human plasma, serum or other biological samples.

Assay performance & characteristics

Human galectin-3 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.082	0
0.625	0.111	0.029
1.25	0.15	0.068
2.5	0.213	0.131
5	0.404	0.357
10	0.782	0.70
20	1.542	1.46
40	2.663	2.581

A. Typical standard curve



Human galectin-3 standard curve (4 parameters)

- B. Sensitivity: 0.145 ng/mL.
- C. Precision: Intra-assay C.V. <5.5%; Inter-assay C.V. < 6.2%
- **D. Recovery:** 74-91%

References:

- 1. Dumic, J. et al. (2006) Biochim Biophys Acta. 1760(4):616-35.
- 2. Wang, X. et al. (2015) Am. J. Alzheimers Dis. Other Demen. 30(8):729-32.
- 3. He, X.W. et al. (2017) Sci. Reports 7: 40994.
- 4. Nangia-Makker, P. et al. (2007) Cancer Res. 67(24):11760-8
- 5. Li, P. et al. (2017) Cell. 167(4):973-984.
- 6. Dong, R. et al. (2018) Int. J. Mol. Med. 41(2):599-614.
- 7. Li, Y. et al. (2020) J Biochem Mol Toxicol. 30: e22463.

28




Human Haptoglobin ELISA Kit (Cat. No. 31400)


INTRODUCTION: Human haptoglobin is an acute-phase glycoprotein produced predominantly by liver. It is composed of two α subunits (M.W. 16-23 kDa) and two β subunits (M.W. 35-40 kDa). Haptoglobin can bind to free hemoglobin released from lysed erythrocytes and prevent the formation of free radical superoxide that can be formed by the reaction of oxygen and iron from hemoglobin. It is also known to be involved in immune regulation and anti-inflammation. Elevated amount of haptoglobin is observed during infections and inflammations, obesity, tissue damage etc. Hence, haptoglobin is used as a biomarker to detect acute allograft rejection, proliferative diabetic retinopathy (PDR) and diabetic kidney disease (DKD). Additionally, low haptoglobin levels are mainly observed during hemolytic anemia.

IMD has developed a highly specific-ELISA for accurate measurements of haptoglobin in human plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve





B. Sensitivity: 1.56 ng/mL

- C. Precision: Intra-assay C.V. <5.63%.; Inter-assay C.V <1%.
- **D.** Spiking: The recovery of human haptoglobin spiked at 5 ng/mL is 104.4%.

E. Linearity:

Sample Dilution	Absorbance (450 nm)	Concentration (ng/mL)	Recovery (%)
1:2	0.601	5.843	95.37
1:4	0.269	6.010	98.10
1:8	0.116	6.790	110.82

- 1. Andersen, C.B. et al. (2012) Nature 489:456.
- 2. Shin, A.W. et al. (2014) Am. J. Hematol. 89(4):443-7
- 3. Körmöczi, G.F. *et al.* (2016) Eur. J. Clin. Invest. **36(3)**:202-9.
- 4. Chiellini, C. *et al.* (2004) J. Clin. Endocrinol. Metab. **89(6)**:2678-83.
- 5. Maffei, M. et al. (2016) Endocr. Rev. 37(4):403-16
- 6. Yang, J.K. et al. (2017) Diabetes Care 40(2):253-260




Human PAI-1 ELISA Kit (Cat. No. 31070)


INTRODUCTION: Plasminogen activator inhibitor-1 (PAI-1), also known as Serpin E1, belongs to serine protease inhibitor family and its primary function is to inhibit tissue-type as well as urokinase-type plasminogen activator. PAI-1 is a secretory protein abundantly expressed in adipose tissue (therefore classified as adipokines), with a molecular weight of 43 kDa. Active form of PAI-1 is unstable and hence undergoes conformational change to inactive form. Circulating PAI-1 level is closely associated with obesity-related insulin resistance, metabolic syndrome, breast cancer, depression and asthma. Additionally, absence of PAI-1 is linked to abnormal bleeding.

IMD has developed monoclonal antibody-based ELISA for specific measurements of PAI-1 in human plasma, serum or other biological samples.

Assay performance & characteristics



A. Typical standard curve



Human PAI-1 standard curve (4-parameter)

B. Sensitivity: 0.031 ng/mL.

C. Precision: Intra-assay C.V. < 3.02%; Inter-assay C.V. < 5.63%.

D. Recovery: 88-114%

E. Specificity:

The antibodies used in this assay are specific to human PAI-1 and do not cross-react with mouse and rat PAI-1, and other cytokines or hormone molecules including human resistin, $TNF\alpha$, ANGPTL4, insulin, leptin and IL6.

- 1. Vousden, K.A. et al. (2019) Sci. Rep. 9:1605
- 2. Vaughan, D.E. et al. (2017) Arteroscler. Thromb. Vasc. Biol. 37:8
- 3. Wong, I.K. et al. (2018) Arterioscler. Thromb. Vasc. Biol. 38(10):2519-2527
- 4. Vague, P. *et al.* (1986) Metabolism **35**: 250-253
- 5. Sobel, Be. et al. (2003) Arterioscler. Thromb. Vasc. Biol. 23: 1979-198.
- 6. Heilbronn, L.K. et al. (2013) PLoS One 8(10): e78864.
- 7. Chen, B. et al. (2006) Biochem. Biophys. Res. Commun. 341(2):549-56



Human Thrombospondin-2 ELISA Kit (Cat. No. 31101)


Introduction: Thrombospondin-2 (TSP-2/THBS-2) is a 150 kDa secreted calcium-binding glycoprotein involved in regulation of proliferation, aggregation, motility, angiogenesis, tumor progression, and wound healing. It has been implicated in regulating cardiovascular inflammation and the immune response and maintaining the integrity and function of cardiac structures. TSP-2 is known to modulate cellular interactions with extracellular matrix, and as such, is considered a matricellular protein. Elevated TSP-2 expression has been observed in many types of cancers. Serum TSP-2 is a promising diagnostic biomarker for early non-small-cell lung cancer and liver fibrosis and is also a predictor for response to treatment with intravenous immunoglobulin in children with Kawasaki disease.

IMD has developed a monoclonal antibody-based ELISA for quantitative determination of TSP-2 in human plasma, serum or other biological samples.

Assay performance & characteristics

Human TSP-2 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.094	0
0.156	0.133	0.039
0.312	0.187	0.093
0.625	0.293	0.199
1.25	0.474	0.38
2.5	0.812	0.718
5	1.317	1.223
10	2.096	2.002

A. Typical standard curve



B. Sensitivity: 0.156 ng/mL.

C. Precision: Intra-assay C.V. < 4.6%.; Inter-assay C.V. < 7.2%.

- 1. Streit, M. et al. (1999) Proc. Natl. Acad. Sci. U S A 96(26): 14888-14893
- 2. Oshika, Y. et al. (1998) Clin. Cancer Res. 4:1785-1788
- 3. LaBell, T.L. & Byers P.H. (1993) Genomics 17:225
- 4. Reinecke, H. et al. (2013) Cardiovasc. Pathol. 22(1):91-95
- 5. Bae, O.N. et al. (2013) Arterioscler. Thromb. Vasc. Biol. 33(8):1920-7
- 6. Abu El-Asrar, A.M. et al. (2013) Acta Ophthalmol. 91(3): e169-77
- 7. Tsai, E.A. et al. (2016) Cell. Mol. Gastroenterol. Hepatol. 2(5):663-675
- 8. Po-Chun, C. et al. (2017) J. Hematol. Oncol. 10:33
- 9. Liu, J.F. et al. (2018) Biochem. Pharmacol. 155:537-546



Human TIMP-1 ELISA Kit (Cat. No. 31103)

INTRODUCTION: Tissue inhibitor of metalloproteinases-1(TIMP-1) is an inducible protein containing 184 amino acids and widely synthesized in many cells. Production of TIMP-1 is induced by cytokines and phorbol esters. It inhibits the activity of matrix metalloproteinases (MMP) which plays a major role in the development of liver fibrosis. The catalytic site of MMP is blocked by the non-covalent binding of TIMP-1. Sites present in the N-terminal domain of TIMP-1 binds to the substrate binding site of MMP and inhibits its function. A higher liver fibrosis rate was observed in TIMP-1 knock out mice after carbon tetrachloride (CCl₄) injection when compared to wild mice. Due to the higher serum level of TIMP-1 in patients with non-alcoholic fatty liver disease (NAFLD), TIMP-1 was categorized as an independent predictor for fibrosis. Additionally, it was also found to inhibit B-cell apoptosis and promote erythropoiesis.

IMD has developed monoclonal antibody-based ELISA for specific measurements of TIMP-1 in human plasma, serum or other biological samples.


Assay performance & characteristics

in Typical Standard Carto				
	Human	Abaarbanaa		

A Typical standard curve

Human TIMP-1 (pg/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.091	0
31.2	0.118	0.027
62.5	0.161	0.07
125	0.257	0.166
250	0.449	0.358
500	0.876	0.785
1000	1.544	1.453
2000	2.494	2.403

Human TIMP-1 standard curve (4-parameter)



B. Sensitivity: 31.2 pg/mL.

C. Precision: Intra-assay C.V.<5.0%; Inter-assay C.V. <3.9%.

- 1. Schultz, R.M. et al. (1988) Cancer Res. 48:5539
- 2. Gomez, D.E. et al. (1997) Eur. J. Cell Biol. 74:111
- 3. Wang, H. *et al*. (2011) Cell Biosci. **1**:14
- 4. Yilmaz, Y. & Eren, F. (2019) Eur. J. Gastroenterol. Hepatol. 31(1):43-46
- 5. Knight, B.E. et al. (2019) Front. Mol. Neurosci. 12:220



Rapid Human Cystatin C ELISA Kit (Cat. No. 31241)


INTRODUCTION: Human Cystatin C (or cystatin 3) is composed of 120 amino acid residues with a molecular weight of ~13.4 kDa. It belongs to the cystatin superfamily, which inactivates lysosomal cysteine proteinases. As a strongly cationic and low-molecular weight (13.4 kDa) protein, it is almost freely filtered across the glomerular membrane, and is therefore used as a biomarker of kidney function. A growing body of evidence suggests that cystatin C is a more reliable biomarker of glomerular filtration rate than creatinine. Additionally, increased serum levels of cystatin C are associated with several types of cardiovascular disease, including myocardial infarction, stroke, heart failure, peripheral arterial disease and metabolic syndrome.

IMD has developed a highly specific antibody-based ELISA for accurate measurements of Cystatin C in human plasma, serum or other biological samples.

Assay performance & characteristics

Cystatin C	Absorbance	Blanked
(ng/mL)	(450 nm)	Absorbance
0	0.069	0
0.23	0.127	0.058
0.46	0.187	0.118
0.93	0.309	0.24
1.87	0.539	0.47
3.75	0.946	0.877
7.5	1.685	1.616
15	2.551	2.482

A. Typical standard curve



Human cystatin C standard curve (4-parameter)

B. Sensitivity: 0.23 ng/mL.

- C. Precision: Intra-assay C.V. 7.9%, inter-assay C.V. 7.2%.
- **D. Recovery:** 94-107%

E. Specificity:

The antibodies used in this assay are specific to human cystatin C and do not cross-react with mouse and rat cystatin C, and other cytokines or hormone molecules.

- 1. Stevens, L.A. *et al.* (2008) Am. J. Kidney Dis. **51**:395-406.
- 2. Zethelius, B. et al. (2008) N. Engl. J. Med. 358 (20): 2107-16.
- 3. Dharnidharka, V.R. et al. (2002) Am. J. Kidney Dis. 40 (2): 221-6.
- 4. Hermida, J. & Tutor, J.C. (2006) Ther. Drug Monit. 28 (3): 326-31.
- 5. Deo, R. et al. (2008) Am. Heart J. 155 (1): 62-8.
- 6. Servais, A. et al. (2008) Am. J. Med. **121 (5)**: 426-32.



Mouse Lipocalin-2 ELISA Kit (Cat. No. 32050)


INTRODUCTION: Refer Cat. No. 31050 for more information.

IMD has developed monoclonal antibody-based ELISA for specific measurements of lipocalin-2 in mouse plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve

Mouse LCN2 (ng/mL)	Absorbance (450 nm)	Blanked Absorbance
0	0.074	0
0.156	0.122	0.048
0.312	0.162	0.088
0.625	0.239	0.165
1.25	0.406	0.332
2.5	0.707	0.633
5	1.236	1.162
10	2.082	2.008



Mouse LCN2 standard curve (4-parameter)

B. Sensitivity: 0.156 ng/mL.

C. Precision: Intra assay C.V. <6.7%; Inter assay C.V. <7.2%

D. Specificity: The antibodies used in this assay are specific to mouse LCN2 and do not cross-react with human LCN2, and other cytokines or hormone molecules.

- 1. Wang, Y. et al. (2007) Clin. Chem. 53(1):34-41
- 2. Law, I.K. et al. (2010) Diabetes 59 (4)
- 3. Xiang, Y. et al. (2011) Diabetes care 34 (7): 1639-41
- 4. Flo, T.H. *et al.* (2004) Nature **432**:917-21.
- 5. Mishra, J. et al. (2005) Lancet 365:1231-6.
- 6. Yang, J. et al. (2009) Proc Natl Acad Sci U S A 106(10) :3913-8.




Human Lipocalin-2 ELISA Kit (Cat. No. 31050)


INTRODUCTION: Lipocalin-2(LCN2), also known as neutrophil gelatinase-associated lipocalin (NGAL), siderocalin or neutrophil lipocalin (NL), is a secretory glycoprotein covalently bound to a matrix metalloproteinase 9 (MMP9) and has a molecular weight of approximately 25 kDa. During acute kidney injury (AKI), NGAL is secreted at high levels into the bloodstream and urine within 2 hours of injury. Therefore, NGAL is a more accurate and sensitive biomarker for diagnosing AKI than serum creatinine. In fact, the increase in urinary excretion of NGAL has been proven to be due to tubular alterations which take place well before any damage that can be detected by other methods. Additionally, lipocalin-2 has been identified as an adipokine closely related to obesity-related insulin resistance, diabetes and fatty liver diseases.

IMD has developed a highly specific monoclonal antibody pair-based ELISA for quantitative measurements of lipocalin-2 in human plasma, serum or other biological samples.

Assay performance & characteristics

A. Typical standard curve





35

B. Sensitivity: 0.39 ng/mL.

C. Precision: Intra-assay C.V. 1.84%; Inter-assay C.V. 6.77%

D. Specificity:

The antibody pair used in this assay is specific to human lipocalin-2 and does not cross-react with mouse and rat lipocalin-2, and other cytokines or hormone molecules.

- 1. Wang, Y. et al. (2007) Clin. Chem. 53(1):34-41
- 2. Law, I.K. et al. (2010) Diabetes 59 (4): 872-82
- 3. Xiang, Y. *et al.* (2011) Diabetes care **34 (7)**: 1639-41
- 4. Flo, T.H. *et al.* (2004) Nature **432**:917-21.
- 5. Mishra, J. *et al.* (2005) Lancet **365**:1231-6.
- 6. Yang, J. et al. (2009) Proc. Natl. Acad. Sci. U S A 106(10):3913-8.



ImmunoDiagnostics Limited

Hong Kong

Unit 513, 5/F, Biotech Centre 2 No.11 Science Park West Avenue, Hong Kong Science Park, Sha Tin, Hong Kong

Email: <u>sales@immunodiagnostics.com.hk</u> Phone: +852 3502 2780

Canada

3330 Bayview Avenue Block #6, Toronto, M2M 3R8 Ontario, Canada

Email: info@immunodiagnostics.ca Phone: +1-437-886-5136

www.immunodiagnostics.com.hk