



## **Monoclonal Mouse Anti-Human IgG, HRP**

**Catalog Number: 221120**

### **Product Details**

Size	1 ml
Species Reactivity	Human
Host/Isotype	Mouse/IgG
Class	Monoclonal
Type	Secondary Antibody
Conjugate	HRP
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Storage Buffer	PBS, pH7.4, 50% glycerol, 0.05%NaN <sub>3</sub>
Storage conditions	Store at 4°C for short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

### **Application**

ELISA	1:4000-1:20000
-------	----------------

### **Product Specific Information**

This antibody can recognize human IgG Fab.