



Monoclonal Antibody against Human CRP (8G1)

Catalog Number: 21120**Size: 100 µg****Host: Mouse**

Introduction to the Molecule

C-reactive protein (CRP) is a circulating protein mainly secreted from the liver. This acute phase protein consists of five identical non-glycosylated subunits of 23 kDa, that give rise to a symmetrically arranged globular protein with molecular weight of approximately 120 kDa¹. It has long been recognized that CRP is closely related to immunology, inflammation and host defence, as a result it has been used as an inflammatory marker. There is accumulating evidence suggesting the important role that CRP plays in mediating cardiovascular diseases (CVD) and type 2 diabetes²⁻⁴. Normally, CRP is present only in a trace amount in circulation (<1 µg/mL)^{5,6} but can increase over 1,000-fold under acute inflammatory state. Individual with blood CRP levels <1 µg/mL, 1-3 µg/mL and >3 µg/mL is considered to have low, moderate and high risk, respectively, of CVD and myocardial infarction⁷. Therefore, blood CRP level has become a promising measure of CVD risk^{8,9}.

Purification

Protein G affinity purification

Immunogen

Recombinant full-length human CRP in *E.coli*.

Specificity

The antibody detects human CRP.

Formulation & Storage

Liquid in phosphate-buffered saline (PBS). Store at -20°C for less than one week. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/thaw cycles.

Application/Usage

This antibody can be used as a capture antibody in a human CRP ELISA in combination with monoclonal anti-human CRP antibody (Cat. No.: 21121) as detection antibody.

References

- [1]Thompson D., Pepys M.B. and Wood S.P. (1999) *Structure*, 7, 169-177.
- [2]Festa A, D'Agostino R. Jr., Tracy R.P. and Haffner S.M. (2002) *Diabetes*, 51, 1131-1137.
- [3]Verma S. and Yeh E.T. (2003) *Am J Physiol*, 285, R1253-R1258.
- [4]Jialal I., Devaraj S. and Venugopal S.K. (2004) *Hypertension*, 44, 6-11.
- [5]Kindmark C.O. (1972) *Scand J Clin Lab Invest*, 29, 407-411.
- [6]Macy E.M., Hayes T.E. and Tracy R.P. (1997) *Clin Chem*, 43, 52-58.
- [7]Ridker P.M. (2004) *Am Heart Hosp J*, 2 (4 Suppl 1), 4-9.
- [8]Benzaquen I.R., Yu H. and Rifai N. (2002) *Crit Rev Clin Lab Sci*, 39, 459-497.
- [9]Pearson T.A. et al., (2003) *Circulation*, 107, 499-511.