



## **Goat anti-Rabbit IgG (H+L) Secondary Antibody, HRP**

**Catalog Number: 221117**


### **Product Details**

Size	1ml
Species Reactivity	Rabbit
Host/Isotype	Goat/IgG
Class	Polyclonal
Type	Secondary Antibody
Conjugate	HRP
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Storage Buffer	PBS, pH7.4, 50% glycerol, 0.05%NaN <sub>3</sub>
Storage conditions	Store at 4°C for short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

### **Application**

Western Blot (WB)	1:4000-1:12000
ELISA	1:4000-1:20000

### **Product Images**



Pancreatic islets and exocrine cells isolated from 10-week-old male C57BL/6N mice were subjected to immunoblotting using an antibody against APPL1 (IMD #11130) or GAPDH (CST #2118).