



## Goat anti-Mouse IgG (H+L) Secondary Antibody, HRP

Catalog Number: 221116


### Product Details

Size	1ml
Species Reactivity	Mouse
Host/Isotype	Goat/IgG
Class	Polyclonal
Type	Secondary Antibody
Conjugate	HRP
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Storage Buffer	PBS, pH7.4, 50% glycerol, 0.05%NaN <sub>3</sub>
Storage conditions	Store at 4°C for short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

### Application

Dot blot	Use at an assay dependent dilution
Western Blot (WB)	1:3000-1:15000
ELISA	1:5000-1:20000

### Product Images



siRNA-transfected INS-1E cells were treated with CKS for indicated time points, followed by immunoblotting using an antibody against APPL1, phospho (p)-IκBα (Ser32/36) (CST #9246), IκBα, p-p65 (Ser536), p65 or GAPDH.