



Mouse Anti-Human IgG (FC), HRP

Catalog Number: 221119

Product Details

Size	1 ml
Species Reactivity	Human
Host/Isotype	Mouse/IgG
Class	Monoclonal
Type	Secondary Antibody
Conjugate	HRP
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Storage Buffer	PBS, pH7.4, 50% glycerol, 0.05%NaN ₃
Storage conditions	Store at 4°C for short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

Application

ELISA	1:5000-1:20000
-------	----------------

Product Specific Information

This antibody can recognize Fc region of human IgG1/2/3/4.